



Name	Walz Phyto-PAM
Description	<p>The PHYTO-PAM can be used to assess the photosynthetic performance and light-adaptation status of the various types of phytoplankton. Furthermore, it allows to determine the content of active chlorophyll in natural surface waters down to 0.1 <math>\mu\text{g Chl/l}</math>, and also to differentiate between differently pigmented groups of algae (green algae, diatoms and cyanobacteria).</p> <p>Fluorescence is excited alternatingly at high repetition rates by <math>\mu\text{sec}</math> pulses of 470, 520, 645 and 665 nm light originating from light emitting diodes (LED). A miniature photomultiplier detector serves for extremely sensitive fluorescence detection. It is protected against light damage by automatic switch-off circuitry.</p> <p>A Spherical Micro Quantum Sensor US-SQS is included to measure the photosynthetically active radiation (PAR) within the two types of cuvettes to calibrate the light list of the PHYTO-PAM.</p> <p>Two PHYTO-PAM system versions with different emitter-detector units are available for different types of samples. They are all based on the same central Power-and-Control-Unit PHYTO-C:</p> <ul style="list-style-type: none"> <li>• MODULAR version, with separate LED-array cones for measuring light and actinic illumination, featuring a 10 x 10 mm cuvette.</li> <li>• COMPACT version as more portable unit for field applications, with all opto-electronical components being mounted in one miniature housing, featuring 15 mm <math>\varnothing</math> round cuvette.</li> </ul>
Services	<p>Measure of photosynthetic performance and light-adaptation status on natural samples</p> <p>Measure of photosynthetic performance and light-adaptation status on microalgal cultures</p> <p>Measure of photosynthetic performance and light-adaptation status in laboratory experiments</p>
Contacts	<p>Maria Saggiomo          Tel. +39 081 5833240          e-mail: maria.saggiomo(at)szn.it</p>

